

Mojia

TITAN

HiFU - RF - EMS

Multifunctional Beauty Device for At-Home Use



In the Box


User Manual can also be downloaded from the Mojia website.



Main Device

Power plug

Getting To Know Your Device



Product Description

Mojia TITAN is a multifunction device, incorporating Hifu, Radiofrequency, and EMS

Hifu - uses ultrasound technology that directly acts on the SMAS layer promoting new fascia cells.

Radiofrequency - can emit 1 million micro-waves per second to penetrate and heat the deeper layer of the dermis, promoting blood flow, collagen production, and helps to restore skin elasticity.

EMS - provides muscle stimulation which helps improve edema, and lift, firm and tone skin.

Blue light indicates: HiFu mode.

Red light indicates: RF mode.

Purple light indicates: EMS mode

Modes Explained



HiFu - High Intensity Focused Ultrasound can channel beams of ultrasound energy precisely up to 4.5 mm bypassing the superficial skin to target the SMAS layer.

This technology significantly tightens and lifts the skin.



RF allows energy penetration up to 3mm to the dermis hereby boosting collagen and elastin which tighten the skin without damaging the top layer.



EMS uses mild current to deeper stimulate the muscle by repeatedly causing it to contract. This significantly promotes muscle movement and improve skin firmness.

How to Use Your Device



- Plug into power outlet. This devices is not battery operated and must be connected to power.
- Press the power button for 3 seconds to turn on.
- The LED screen will display Intensity level, Timer, and Mode.
- The first mode to display is HiFu. Toggle between modes by pressing 'M' button.
- Press the 'L' button to change the intensity level. There are 5 levels which are reflected by a bar indicator at the top left of the screen.
- Hold the power button down until device turns off.
- Disconnect from power and clean with a damp cloth,



WARNING: Ensure power outlet is not switched on whilst plugging the cable into your device.

Once cable is connected to your device you can then switch on mains power.

Plugging the cable into your device whilst mains power is switched on at the power outlet may cause damage to your device.

Usage Techniques

Always perform a treatment on clean skin.

Remove makeup, oils and dirt with an oil free cleanser and pat dry.

Remove jewellery that may come in contact, or interfere with the device.

It is recommended to use a quality Conductive Gel

We recommend to use our Electra Conductive Booster Gel for all modes.

ALWAYS be sure to start on Level 1 intensity when using for the first time.



When doing a full face treatment, divide your face into four areas/sections, as per the diagram.

Apply gel liberally to one area at a time.

Divide neck into two, making sure to avoid the thyroid area.



HIFU - RF - EMS **Gliding Movements**

Hifu - Move very slow, approximately 2 to 3 seconds per centimetre (1/4").

RF & EMS - Glide less slowly.

You may do up to three passes over the same area.

If it feels dry whilst gliding add more gel.



RF **Circular Movements**

Using small circular movements.

CAUTION: DO NOT hold in one position when using RF. This may result in skin damage and burns.



HIFU - EMS **Holding / Stamping Movements**

Use the stamping method by holding the device in one place - referred to as 'shot'. Each shot must be no longer than 3 seconds.

Use gel liberally and reapply often if doing this method.

You may do up to three passes.

DO NOT USE THIS METHOD FOR THE RF MODE



WARNING

Do not use on the upper eyelids, eyeballs, around the ears, or over the thyroid area.

Do not use on mucosa - a membrane that lines various cavities in the body.

Pre-Programmed Duration

MODE	DURATION per treatment	FREQUENCY Week 1 - 6	FREQUENCY Week 6+
HIFU	10 minutes	3 x weekly	1 x weekly
RF	10 minutes	3 x weekly	1 x weekly
EMS	10 minutes	3 x weekly	1 x weekly

IMPORTANT - Treatment frequencies in the above table are based on the recommended duration time of 10 minutes per **treatment. NOT** per **area/section**

It is recommended to use the modes from deeper layers to epidermis. Hifu / RF/ EMS in that order.

You can use one mode for a longer duration, for example: Hifu for 20 minutes. When doing this the remaining modes must be not used for longer than 5 minutes each.

You can use one mode only but this single use mode must not exceed 30 minutes.

Example: Hifu 30 minutes only - NO RF or EMS

You can be flexible with the duration of modes to suit your needs however, you must stay within the recommended total duration of 30 minutes.

Frequency of use must not exceed the recommendations.

Product Specifications

HIFU

Frequency: 3 MHz
Power: 5J/cm²
Active area: 3cm²
Depth: max 4.5mm

RF

Frequency: 2 MHz
Power: 4J/cm²
Active area: 3cm²
Depth: max 3mm

EMS: Pulsed Current

Blue light: represents hifu mode

Red light: represents RF mode

Purple light: represents EMS mode

Power rating: 15w

Power consumption: 500 mAH

Output Frequency: 50-60 Hz

Net Weight of main unit only: 164gm













Contraindications

Do not use if you have any of the following, but not limited to -

- Heart disease, abnormal blood pressure, malignant tumors, cancer
- Migraines
- Seizures
- Pacemaker
- Currently undergoing dental treatment
- Dermatitis, psoriasis, or other skin disease
- Pregnant
- Have had recent chemical peels or similar salon treatments
- Botox or filler within 2 weeks

If unsure, please consult your practitioner for further advice.

Precautions and Warnings

-  Do not use for any other purpose than it's intended use as recommended in this manual.
-  Do not use if the device is defective as this can lead to personal injury.
-  Do not use the adapter with voltage excess of your normal electricity voltage (either 100V or 240V) to avoid damage.
-  Do not use if the power cord is damage.
-  Plug into power outlet and into device before switching power on.
-  After using, turn off the device and switch power off at the power outlet.
-  Remove the plug from the outlet by holding the plug, do not pull the cord to remove as this may damage the cord and cause injury.
-  A certain amount of heat when using this device is to be expected however, if you feel extreme or abnormal heat cease use immediately.
-  Do not disassemble, modify or repair this device. Please contact Mojia for assistance .
-  The device must not come in contact with water during treatment. This may lead to serious injury or electric shock. It may also cause fire or cause the device to malfunction.
-  This device is not waterproof - do not submerge device in water or get it wet.
-  Store the device in its box in a dry place.



Start on Level 1 when using this device for the first time.



Do not press the head of the device hard on the skin.



Do not stay on the same area for longer than 3 seconds, otherwise this may cause skin problems.

T



Do not use around the soft tissue area under the eyes, on the eyelid, or over the thyroid.



If you should feel any abnormal pain when using this device, stop immediately and seek medical advice.



Mojia Australia

www.mojiaaustralia.com

Email: info@mojiaaustralia.com



